TRANSLATION STAGE

CONTROL

Present by Nunnalin Auephattanaphong (Internship)

Robotics and AI Engineering, CMU.

ABOUT INTERNSHIP PROJECT (1)

Get to know instruments

2

TOPICS

Find Solution

3

Experiment & Application

4

Conclusion

5

O1 About Internship project

Translation stage control project is a mini project that part of EVWaCo project from NARIT.

Evanescent Wave Coronagraph or EvWaCo is project that study optical instruments for study about binary star system or exoplanetary system . EvWaCo aim to develop optical instruments for separate light from bright star to observe the companion star


The main idea of translation stage control project is control translation stage moving light source along horizontal axis to make the light source closest align to imaging system for get high quality image by using PID closed loop control.

Why we use PID Control?

- Stability and Responsiveness : The D term dampens oscillations, while the I term corrects steady–state errors, leading to a more stable system.
- Improved Accuracy : The I term ensures the system output matches the desired setpoint more precisely.
- Flexibility: PID provides three tunable parameters (Kp, Ki, Kd) to optimize performance for different systems and dynamic behaviors.

In summary, PID control offers balanced error management, enhancing stability, accuracy, and responsiveness over pure error control.

- Stands for "Proportional, Integral, Derivative".
- It is closed-loop feedback control.
- instrument used by control engineers to regulate temperature, flow, pressure, speed, and other process variables in industrial control systems.



https://microcontrollerslab.com/pid-controller-implementation-using-arduino/

P: Proportional

- tuning -> correcting a target proportional to the difference
- based on the current error (setpoint measured output of the system)

I: Integral


- tuning -> cumulating the error result from the "P" action to increase the correction factor
- <u>Increase</u> -> increases the <u>contribution of the accumulated error</u> over time to the control signal

D: Derivative

- tuning -> minimize overshoot by slowing the correction factor
- reduce the rate of change of the error
- helps improve the stability and responsiveness of the control loop


P Controller

- Use only error term to affect the output
- simply multiply e(t) * Kp
- provides a linear relationship between the error of a system and the controller output of the system




PI Controller

- a combination of Proportional controller action and Integral controller action
- correlates the controller output to the error and the integral of the error.
 provides a faster response time than I-only control
- Eliminates steady state error




PD Controller

- combination of feedforward and feedback control
- contains proportional control's damping of the fluctuation
 - derivative control's prediction of process error output
- Improves transient response, reduces overshoot



PID Controller

- a combination of all three types of control methods
- most commonly used
- balances steady state and transient response



O2 Get to know instruments

INSTRUMENTS

Mainly focus on

- KDC101 -> Controller for control DDS100/M
- DDS100/m -> One axis translation Stage
- ASI1600mm PRO -> Camera for capture images



chrome- KDC101

extension://efaidnbmnnnibpcajpcglclefindmkaj/https://optome chs.com/pdf/manuals/KDC101/KDC101_ManualforKinesis.pdf



DDS100/M

 $\label{lem:https://www.thorlabs.com/thorproduct.cfm?partnumber=DDS100/M$




ASI1600mm PRO

https://telescopes.net/zwo-asi1600mmpro-monochrome-cmos-camera.html

Optical components


- Light source
- Mirrors
- lens

1. Find error of translation stage




- Moving translation stage into steps for find error of stage .
- Home position can't change.
- Home position is 0.0000 .
- Error of translation stage is <u>0.0100 0.0500</u>.

2. Capture image from camera




SETUP



O3 Find Solution

WORKFLOW



mainly focus on


- Find distance between 2 circles
- Find solution of overlapped circles condition

1. Find center coordinates of circle by

- Using OpenCV for find edge , contour area and coordinate of center point.
- Use Matplotlib for plot circle.


2. Find distancebetween 2 circle by

- 1.Find coordinate of center of 2 circles
- 2.Subtrac coordinate of center of 2 circles




3. Find center coordinates and distance of circle

METHOD 1



Can't use this condition with image that has 2 circles that are overlapping

3. Find center coordinates and distance of circle




- Coding program for draw and plot circles from circle formula .
- Initialize as 2D array .
- Plot histogram for see intensity of circle

This code works for overlapped circle condition.

METHOD 3

3. Find center coordinates and distance of circle




main idea of this concept is


Use bitwise operation (AND, XOR) for find contour area that not overlapping. Then draw bounding box around the area then find height and width then calculate center of the area from height and width. Then find distance between 2 center point.

METHOD 3

3. Find center coordinates and distance of circle



OVERVIEW MAP




METHODS	ADVENTAGES	DISVENTAGES
METHOD 1	Easy for operate	Can't use this method with overlapped images
METHOD 2	Can use this method with overlapped images	Can't find distance between 2 images.
METHOD 3	Can use this method with overlapped images	Take more time to process


```
center of ref - x : 2253.0 , y : 1672.5
center of object - x : 2341.0 , y : 1861.5
Distance between objects - x : -88.0 , y : -189.0
-----Normalize-----
CX ref = 8.5614
CY ref = 6.3555
center x = 8.8958
center Y = 7.0737
disX = -0.3344
disY = -0.7182
```

Normalized data from pixel to mm.


04 Experiment & Application

EXPERIMENT





Test code in lab with method 3 as solution for collect real data and find error of process of code when run in real situation.







PROBLEM: Have only one light source so it can't see 2 dot in one image at the same time.

Solution


- Find centroid and position of dot both in reference image and collected image.
- Compare position of current image with reference image.

(Still use bitwise operation method in code for overlapped condition)





EXPERIMENT





<u>PROBLEM</u>: Program doesn't read current images that had capture. It's read and overwrite old images in directory instead.

Solution: Coding new program that has infinity loop and can read images from current images.

```
for path in Dir Read('s'):
   print('Capturing image')
   if i < 10:
       filename = '00'+ str(i)+' image lab.png'
       filename = '0'+ str(i)+' image lab.png'
   camera.set image type(asi.ASI IMG RAW16)
   camera.capture(filename=save path+filename)
   print('Saved to %s' % filename)
   print("----")
   time.sleep(0.5)
   disX = Draw Contour(path)
   err = PID(Decimal(1) , Decimal(0.08), Decimal(0.01) , reference , Decimal(disX))
   print("Error : " + str(err))
   if err > reference :
       new position = pos-err
       kcube.MoveTo(new_position, 7000)
      print("New position : " + str(new position) )
   elif err < reference:
       new position = pos+err
       kcube.MoveTo(new position, 7000)
       print("New position : " + str(new position) )
   elif err == reference:
       break
   time.sleep(0.5)
   error.append(disX)
   plt.plot(error)
   plt.gca().invert yaxis()
   plt.show()
   i = i+1
```




PROBLEM: Program read and draw wrong contour area. It's read the random tiny pixels instead of the real reference pixels. So, it's effect to calculate distance between current pixels and reference pixels.

Solution: Add max(contours1, key=cv2.contourArea) to find and draw the biggest contour as reference pixels

```
if len(contours1) > 0 and len(contours2) > 0:
    cnt1 = max(contours1, key=cv2.contourArea)
    cnt2 = contours2[0]

    x_ref,y_ref,w_ref,h_ref = cv2.boundingRect(cnt1)
    x,y,w,h = cv2.boundingRect(cnt2)
```

Application have 2 versions



```
DetaultValue: 20
Description: 'Sensor temperature(degrees Celsius)'
IsAutoSupported: False
IsWritable: False
MaxValue: 1000
MinValue: -500
Name: 'Temperature'

Do you want to change paramter ? (Y/N) : Y

Please enter kp value : 34

Please enter ki value : 2

Please enter kd value : 0.2

Please enter first position : 70
```


UI Version

Command line Version

Command line version for users.

```
def main():
              answer = input("Do you want to change paramter ? (Y/N) : ")
              if answer == "Y" :
                  kp = float(input("Please enter kp value : "))
                  ki = float(input("Please enter ki value : "))
                  kd = float(input("Please enter kd value : "))
                  pos = float(input("Please enter first position : "))
                  new position = pos
                  new pos = []
145
              elif answer == "N" :
                  kp = 35
                  kd = 2.5
                  ki = 0.1
                  pos = 60
                  new position = pos
                  new pos = []
                  print("Please enter the answer (Y/N)")
              for path in Dir_Read('s', path=save_path):
                  print('Capturing image')
PROBLEMS OUTPUT DEBUG CONSOLE TERMINAL PORTS
exe c:/Users/Asus/Desktop/INTERNSHIP PROGRAM/test ui.py
PS C:\Users\Asus\Desktop\INTERNSHIP PROGRAM> & C:\Users\Asus\AppData/Local/Programs/Python/Python310/python.exe c
PS C:\Users\Asus\Desktop\INTERNSHIP PROGRAM> & C:\Users\Asus\AppData/Local/Programs/Python/Python310/python.exe
ead img pid.py"
Do you want to change paramter ? (Y/N) : Y
Please enter kp value :
```


• First version of UI version of application



Have 2 options .


Option 1: Enter values with default values.

• Click "DEFAULT" button then program will start immediately.




Option 2: Enter values with other values.

• Enter each value in each entry box the click "ENTER" button then click "START" for starting program.





Collected feedback from user (researcher)




- Add the title of application
- Display properly dynamic range of PID parameters
- Create Stop and Home button
- Initialize the PID parameters in entry boxes
- Short describe for effect of each PID parameters

Current Version of application for users






Add "Warning Window" for warning when put some values that out of range.



05 Conclusion

Both UI version and input command in terminal version work successfully.

Can control translation stage with PID closed loop control and collect images.



DetaultValue: 20

Description: 'Sensor temperature(degrees Celsius)'
IsAutoSupported: False
IsWritable: False
MaxValue: 1000
MinValue: -500
Name: 'Temperature'

Do you want to change paramter ? (Y/N) : Y
Please enter kp value : 34
Please enter ki value : 2
Please enter kd value : 0.2
Please enter first position : 70


UI Version


Command line Version

Default values of PID are: P: Proportional: 10

I: Integral : 0.1
D: Derivative : 0.2
Initial position : 60

Application can run with the same values and the results still the same. The result is the position of current image is in the same place of reference image with 0.00 error even though change the reference position .





Reference position = 53.37

Reference position = 55.48